


Getting Started

5

The PSR Tool Kit is now in Tablet Mode

This is the desired mode to use the PSR Tool kit in when enrolling patients or doing any of the form and report related administration tasks.



6

Using the Pen & Tablet

The tablet is used with the special graphic pen. You may use the pen for signatures and to fill in form fields. In order for the tablet to recognize form input you must write in Capital Letters.

The tablet is smart and will fill in much of the information needed for your tasks automatically. You may also use the pen to type on a Virtual Keyboard like the one that appears for you at Sign In.

Locate the **Release** button near the HP logo on the frame surrounding the tablets screen.

A. Click the Release button

B. The pen is now partially ejected & can be removed from it's slot just under the HP logo.

