

Soft Tissue Management

Implant Level Impression

Healing abutment

S	NC			RC		
	conical	bottle shape		conical	bottle shape	
D	3.6 mm	4.8 mm	3.3 mm	4.5 mm	6.0 mm	4.4/4.7 mm
H	2.0 mm	3.5 mm	5.0 mm	2.0 mm	4.0 mm	6.0 mm

Customizable healing abutment

	NC	RC
	024.2270	024.4270

Temporary abutment

	NC	RC
	024.2370	024.4370

Temporary abutment

	NC	RC
Crown	024.2371	024.4371
Bridge	024.2375	024.4375

Impression

Impression post

	NC	RC
closed tray	025.2201	025.4201
open tray, short	025.2202	025.4202
open tray, long, 24.0 mm	025.2205	025.4205

Implant analog

	NC	RC
	025.2101	025.4101

Bite registration aid

	NC	RC
8.0 mm	025.2208-04	025.4208-04
12.0 mm	025.2212-04	025.4212-04

LOCATOR® Impression Coping

	NC	RC
	048.197V4	

LOCATOR® analog

	NC	RC
	025.0120-04	

Planning

Anatomic

PLAN abutment

GH	NC		RC	
	0°	15°	0°	15°
2.0 mm	025.2601-04	025.2603-04	025.4601-04	025.4603-04
3.5 mm	025.2602-04	025.2604-04	025.4602-04	025.4604-04

Meso

Abutment*

	NC	RC
	022.2202	022.4202

Gold

PLAN abutment

	NC	RC
Crown	025.2610-04	025.4610-04
Bridge	-	025.4630-04

Abutment*

	NC	RC
Crown	022.2410	022.4410
Bridge	-	022.4430

Customized

Straumann CAD/CAM Abutment*

	NC	RC
Ti and ZrO ₂	order at: http://cares.straumann.com	

Bar

Abutment*

	NC	RC
Gold	022.2701	022.4701
Ti	022.2702	022.4702

Final Restoration

LOCATOR®

PLAN abutment

GH	NC		RC	
	-	-	-	-
1.0 mm	-	-	-	-
2.0 mm	025.2642-04	025.4642-04	025.2642-04	025.4642-04
3.0 mm	025.2643-04	025.4643-04	025.2643-04	025.4643-04
4.0 mm	025.2644-04	025.4644-04	025.2644-04	025.4644-04
5.0 mm	025.2645-04	025.4645-04	025.2645-04	025.4645-04
6.0 mm	025.2646-04	025.4646-04	025.2646-04	025.4646-04

Abutment

	NC	RC
1.0 mm	-	022.4501
2.0 mm	022.2502	022.4502
3.0 mm	022.2503	022.4503
4.0 mm	022.2504	022.4504
5.0 mm	022.2505	022.4505
6.0 mm	022.2506	022.4506

Cementable

PLAN abutment

GH	NC		RC	
	3.5 mm	5.0 mm	5.0 mm	6.5 mm
4.0 mm	1.0 mm	2.0 mm	2.0 mm	3.0 mm
AH	1.0 mm	2.0 mm	2.0 mm	3.0 mm
	2.0 mm	3.0 mm	3.0 mm	4.0 mm

Abutment*

D	NC		RC	
	3.5 mm	5.0 mm	5.0 mm	6.5 mm
4.0 mm	code █	1.0 mm	2.0 mm	3.0 mm
AH	code █	2.0 mm	3.0 mm	4.0 mm
	code █	3.0 mm	4.0 mm	5.0 mm

Multi-Base

PLAN abutment

GH	NC		RC	
	0°	25°	0°	25°
1.0 mm	025.2651-04	025.4651-04	025.2651-04	025.4651-04
2.5 mm	025.2652-04	025.4652-04	025.2652-04	025.4652-04
4.0 mm	025.2654-04	025.4654-04	025.2654-04	025.4654-04

Abutment*

	NC	RC
Crown	022.2731	022.4731
Bridge	022.2732	022.4732

Copings and Auxiliaries

Basal screw (replacement)

	NC	RC
	025.2900	025.4900

Basal screw (replacement)

	NC	RC
	025.2900	025.4900

Basal screw (replacement)

	NC	RC
Crown	025.2900	025.4900
Bridge	-	025.4902

Auxiliaries (replacement), for customized abutment

	NC	RC
Wax up sleeve, for single use	025.2903	025.4903
Basal screw (for titanium)	025.2900	025.4900
Basal screw (for ceramic)	025.2906	025.4906

Basal screw (replacement)

	NC	RC
	025.2900	025.4900

LOCATOR®

Auxiliaries, for LOCATOR® abutment

	NC	RC
Male processing package	048.189V2	
Replacement male	048.190V4-048.194V4	

Coping

D	NC		RC	
	3.5 mm	5.0 mm	5.0 mm	6.5 mm
Crown	023.2316	023.2326	023.4326	023.4336
Bridge	023.2310	023.2320	023.4320	023.4330

Basal screw (replacement)

	NC	RC
	025.2908	025.4908

Transfer aid

D	NC		RC	
	3.5 mm	5.0 mm	5.0 mm	6.5 mm
	025.2911-04	025.2912-04	025.4912-04	025.4913-04

Abutment Level Impression

Impression taking on implant level possible (see left)

Protective cap

D	NC		RC	
	3.5 mm	5.0 mm	5.0 mm	6.5 mm
4.0 mm	024.2312-04	024.2322-04	024.4322-04	024.4332-04
5.5 mm	024.2316-04	024.2326-04	024.4326-04	024.4336-04

Temporary coping

AH	NC		RC	
	4.0 mm	5.5 mm	5.0 mm	6.5 mm
	024.2346	024.2356	024.4356	024.4366

Impression cap

D	NC		RC	
	3.5 mm	5.0 mm	5.0 mm	6.5 mm
4.0 mm	code █	025.2203-04	025.2204-04	025.4203-04
5.5 mm	code █	-	025.2220-04	025.4230-04

Implant analog

AH	NC		RC	
	4.0 mm	5.5 mm	5.0 mm	6.5 mm
	code █	025.2112	025.2122	025.4132
	code █	025.2116	025.2126	025.4136

Abutment Level Impression

Impression taking on implant level possible (see left)

Protective cap**

D	NC			RC		
	0°	25°	0°	25°	0°	25°
3.5 mm	024.2331-04	024.2341-04	024.2306-04	024.4341-04	024.4361-04	024.4306-04

Temporary coping**

D	NC			RC		
	3.5 mm	4.5 mm	4.0 mm	4.5 mm	6.5 mm	4.0 mm
	024.2332	024.2342	024.2307	024.4342	024.4362	024.4307

Impression cap/post

D	NC		RC	
	0°	25°	0°	25°
3.5 mm	025.2232	025.2242	025.2206	025.4242
Post	025.2232	025.2242	025.2206	025.4242
Cap	025.2231-04	025.2241-04	-	025.4241-04

Implant analog

D	NC		RC	
	0°	25°	0°	25°
3.5 mm	025.2131	025.2141	025.2106	025.4141
4.0 mm	025.2131	025.2141	025.2106	025.4141
4.5 mm	025.2131	025.2141	025.2106	025.4141

Abutment Level Impression

Impression taking on implant level possible (see left)

Protective cap**

D	NC			RC		
	0°	25°	0°	25°	0°	25°
3.5 mm	024.2331-04	024.2341-04	024.2306-04	024.4341-04	024.4361-04	024.4306-04

Temporary coping**

D	NC			RC		
	3.5 mm	4.5 mm	4.0 mm	4.5 mm	6.5 mm	4.0 mm
	024.2332	024.2342	024.2307	024.4342	024.4362	024.4307

Impression cap/post

D	NC		RC	
	0°	25°	0°	25°
3.5 mm	025.2232	025.2242	025.2206	025.4242
Post	025.2232	025.2242	025.2206	025.4242
Cap	025.2231-04	025.2241-04	-	025.4241-04

Implant analog

D	NC		RC	
	0°	25°	0°	25°
3.5 mm	025.2131	025.2141	025.2106	025.4141
4.0 mm	025.2131	025.2141	025.2106	025.4141
4.5 mm	025.2131	025.2141	025.2106	025.4141

Abutment Level Impression

Impression taking on implant level possible (see left)

Protective cap**

D	NC			RC		
	0°	25°	0°	25°	0°	25°
3.5 mm	024.2331-04	024.2341-04	024.2306-04	024.4341-04	024.4361-04	024.4306-04

Temporary coping**

D	NC			RC		
	3.5 mm	4.5 mm	4.0 mm	4.5 mm	6.5 mm	4.0 mm
	024.2332	024.2342	024.2307	024.4342	024.4362	024.4307

Impression cap/post

D	NC		RC	
	0°	25°	0°	25°
3.5 mm	025.2232	025.2242	025.2206	025.4242
Post	025.2232	025.2242	025.2206	025.4242
Cap	025.2231-04	025.2241-04	-	025.4241-04

Implant analog

D	NC		RC	
	0°	25°	0°	25°
3.5 mm	025.2131	025.2141	025.2106	025.4141
4.0 mm	025.2131	025.2141	025.2106	025.4141
4.5 mm	025.2131	025.2141	025.2106	025.4141

☐ = CrossFit® Connection NC = Narrow CrossFit RC = Regular CrossFit S = Shape D = Diameter ↙ = Angulation █ = MO 0 (white) █ = MO 1 (shaded)
L = Length H = Height GH = Gingiva Height AH = Abutment Height *Basal Screw included **Occlusal Screw included *04 = Pack of 4

IPS e.max® is a registered trademark of Ivoclar Vivadent, AG. LOCATOR® is a registered trademark of Zest Anchors, Inc. Straumann® and CrossFit® Connection are trademarks or registered trademarks of Straumann Holding AG or its affiliates. All rights reserved.

STRAUMANN® BONE LEVEL IMPLANT PROSTHETIC OPTIONS CHART



Temporary Abutment
Individual esthetics easy to achieve

Gold Abutment for crown
Versatile solution for a wide range of applications

Gold Abutment for bridge

Straumann® CADCAM Abutments
Straumann customized implant prosthetics, ZrO₂ and titanium



Straumann® Anatomic IPS e.max® Abutments
Esthetics you dream of – predictability you deserve, 0° and 15° angled

Anatomic Abutment
Individual esthetics easy to achieve




Cementable Abutment
Prosthetic handling closest to a natural tooth prep

Abutment for Bars
Userfriendly solution for bar-retained dentures



Multi-Base Abutment
Flexible base for screw-retained multiple-unit restorations

LOCATOR® Abutment
Flexible fixation of full dentures



www.straumann.com

<p>International Headquarters Institut Straumann AG Peter Merian-Weg 12 CH-4002 Basel, Switzerland Phone +41 (0)61 965 11 11 Fax +41 (0)61 965 11 01</p>	<p>Straumann USA Straumann USA, LLC 60 Minuteman Road Andover, MA 01810 Phone 800/448 8168 978/747 2500 Fax 978/747 2490 www.straumannusa.com</p>	<p>Straumann Canada Straumann Canada Limited 3115 Harvester Road, 1st Floor Burlington, ON L7N 3N8 Phone 800/363 4024 905/319 2900 Fax 905/319 2911 www.straumann.ca</p>
---	--	---

Straumann products are CE marked 2/10 USIT 231 © Straumann USA, LLC 2008. All rights reserved. Ivoclar Vivadent® and IPS e.max® are registered trademarks of Ivoclar Vivadent, AG.

COMMITTED TO
SIMPLY DOING MORE
FOR DENTAL PROFESSIONALS™

IPS e.max® is a registered trademark of Ivoclar Vivadent, AG.
LOCATOR® is a registered trademark of Zest Anchors, Inc.

IPS e.max® is a registered trademark of Ivoclar Vivadent, AG.
LOCATOR® is a registered trademark of Zest Anchors, Inc.