

LAUREN PUTTING ONE CONTACT IN EYE – 44 SEC.

55-59 s	4	takes right lid off case
59-1.06 s	7	adjusts right contact on her finger
1.06-1.08s	2	picks at contact
1.08-1.19	11	puts lens back into case and swishes around
1.19-1.27	8	tries several times to get contact out of case
1.27-1.31	4	adjusts right contact on her finger
1.32-1.38	6	puts right contact in eye

ISSUES –

Takes long time to situate contact on her finger
Doesn't see dirt on contact and has to clean off again
Takes 8 seconds just to get the contact out of the case

LAUREN TAKING CONTACTS OUT OF EYE ROUND 2

-2-2s	4	unscrews lids of case
2-10s	8	rinses off case in water
10-27	17	opened solution, filled the case with it
27-34	7	took right contact out of eye
34-36	2	put contact into case
36-43	7	took left contact out of eye
43-46	3	puts contact into case

ISSUES –

Rinses case with water and uses immediately
Takes a very long time to fill case with solution